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ENGLISH

All implants Axelmed Paradigma are compatible with all instruments, prosthetic and impression components from code 82001 to code 85050 (see current catalog). All implants Axelmed Paradigma are made of titanium Gr. 4, in accordance with the requirements ISO 5832-2. The Axelmed Paradigma line consists of dental implants, with an internal hexagonal connection and a conical body. The surface is machined at the level of the implant neck, while the remaining part is treated with a blasting process using aluminum oxide and then a triple acidification process. After the surface modifications implants are washed in a ultrasonic bath, treated with dry steam vaporization at high pressure, decontaminated with double cycle of cold Argon plasma, packaged in Class A Bio-Hazard Hood and sterilized by irradiation with Gamma Rays at minimum 25kGy. The prosthetic components are made of titanium (Gr.5) according to ISO 5832-2. Calcinable abutments, impression transfers and accessories (label: Plastic) are made of PMMA. Qualified doctor can perform both submerged surgical technique (two surgical phases) and non-submerged (one-stage surgery). The Axelmed Paradigma implant line can be used only by qualified professionals and it is designed as part of a unique system, the Axelmed Paradigma System: implants, prosthetic components and surgical instruments must be original and used in accordance with the precise recommendations for uses supplied by the manufacturer. The detailed information on the material, the LCS (washing processes, packaging, sterilization), the surgical and prosthetic steps are available in the "Surgical and Prosthetic Axelmed Paradigma Guide". Recommendations: The surgical procedure for the insertion of Axelmed Paradigma implants is complex and highly specialized. So the operator must have completed appropriate training courses. Each implant site should undergo an appropriate diagnostic evaluation, both clinical and radiological. Incorrect procedures can cause loss of the implant and biological damage. The Axelmed Paradigma must be used with instruments specifically designed for oral implantology and with specific prosthetic components of the Axelmed Paradigma line. Pain on percussion, bone infection and bone loss are indicators of implant failure. Packaging: The Axelmed Paradigma implant is packaged in an "IMPLANT KIT", which includes: Alloplastic conical implant made of titanium with surface modification • Titanium box • Titanium Closure cap • Titanium healing cap • First PVC tube containing the implant • First PVC tube cap • Second PVC tube containing the first PVC tube • Cap for the second PVC tube • Label to be placed on the patient's medical record • Packaging box • Package insert inside the box for the shipping. The Axelmed Paradigma implants is packaged also in "ALL-IN-ONE KIT", which includes: Alloplastic conical implant made of titanium with surface modification • Titanium box • Titanium Closure cap • Titanium healing cap • First PVC tube containing the implant • First PVC tube cap • Second PVC tube containing the first PVC tube • Cap for the second PVC tube • Label to be placed on the patient's medical record • 00 Abutment with titanium screw (gold-colour) • Plastic direct transfer • Titanium laboratory screw (blue-colour) • Castable abutment • Laboratory analog • PVC tubes for prosthetic components • Package insert inside the box for the shipping. Indications: The Axelmed Paradigma implants

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Rev.2 del 09.11.2019

can be used by qualified doctors or dentists in order to replace, according to the usual surgical techniques, missing dental roots in the maxillary and mandibular bones. The final aim is the aesthetic-functional rehabilitation of the mouth. The use of Axelmed Paradigma allows to preserve the natural elements adjacent to the missing tooth and guarantees a valid aesthetic and functional solution. This solution allows to realize fixed prostheses in patients with no anchor points "bridge" (therefore forced to use removable dentures). Endosseous dental implants, compared with mobile solutions, permit greater functional stability and bone atrophy reduction. Patient Selection: The prosthetic long-term success of endosseous implants is obtained by a proper integration into the bone, a precise anatomical and surgical planning and, of course, from an evaluation of relative and absolute indications that the dentist should know before the surgery procedure, such as: Psychological assessment of patients; Possible alternatives to implant surgery; General medical history: general health of the patient, metabolic and cardiovascular disease, alterations in intercellular processes and / or alteration of the immune system that can be relative or absolute contraindications for surgery procedure; Local history and healthy of hard and soft tissue; Morphological characteristics (bone quality and quantity); Motivation and maintaining oral hygiene; The overall design of the two dental arches: elimination of microbial contamination and consideration of the masticatory aspect; Specific radiological exams: intraoral RX, RX orthopantomography and computed tomography; Dental models analysis; Absolute Contraindications: Recent myocardial infarction <6 months: Surgical Stress => Tachyarrhythmia - Hypertension - Vasoconstriction controlled - Infection risk for the myocardium, endocardium and pericardium; Severe kidney disorders: impossibility for tubules to reabsorb Ca²⁺ with metabolic loss - Malfunction of PTH - Excess of plasma endotoxin - Risk of infection - Metabolic osteopenia; Osseomalacia: risk of failure for the implant integration; Diabetes treatment-resistant (uncompensated): Macro / Micro systemic angiopathy - Metabolic disorders (Cl, Na, Mg, etc.) - Ready to tissue degeneration - Healing difficulty - Increased risk of infections; Recent radiation therapy with high doses: Reduction of defense mechanisms - Impairment of vascular defenses - Reduction of osteoinduction - Reduction of osteoconduction - Tissue necrosis - Loss of physiological activity of the periosteum - Increased risk of infections; Coagulation disorders: micro-macro bleeding; Advanced liver failure: Hypoproteinemia - Difficulty healing; Chronic alcoholism: liver disorders, cirrhosis - bone aplasia, coagulation disorders, myocardial stress, aneurysms, etc. - Delayed healing - Psychological disorders - Poor hygiene - Increased risk of infections; Drugs: Loss of the sense of priorities - Undernutrition - Predisposition to infections - Psychological disorders - Poor hygiene; Disease and cancer in terminal phase; Relative Contraindications: Chemotherapy; Mild kidney disorders; Liver disorders; Endocrine disorders; Psychological disorders and psychosis; AIDS, HIV-positive cases; Prolonged use of corticosteroids; Disorders of calcium-phosphorus metabolism; Preoperative Measures: The dental office, surgery and operating rooms must be properly equipped with instruments and drugs as required by law. It is necessary to perform the surgery procedure under sterile conditions, thoroughly disinfecting the oral mucosa, oral and peri-oral area and isolating the operative field. For anxious patients, it is recommended the administration of benzodiazepines at least one hour before surgery. And it is advisable to accompany the patient immediately after his discharge. Special Preoperative Measures: Patients at risk: Based on the clinical history and the physical examination, the patient should be classified according to the metabolic or cardiovascular risk. Particularly, perioperative surveillance should be reserved in some patients at high risk due to the following disease: 1) Cardiac or coronary insufficiency 2) Hypertension 3) Diabetes 4) Epilepsy 5) Allergic diathesis 6) Coagulopathies. The key vital signs (E.C.G., P.A., pulse oximetry) should be monitored. The sedoanalgesia is normally recommended for the patient's comfort (dental phobia and leipothymics) and operator surgery. Surgery Procedure: The surgical steps for the use of Axelmed Paradigma should be carefully respected following the protocols illustrated on "Surgical and Prosthetic Axelmed Paradigma Guide". It is advisable to irrigate the surgical field with sterile physiological solution in order to eliminate any residue below the mucoperiosteal flap. Immediate Postoperative Measures: The patient should thoroughly disinfect the oral cavity with suitable solutions. Following surgery, it is advisable to administer to the patient an anti-inflammatory drug. After checking very carefully the effectiveness of hemostasis, the patient can be discharged in accordance with the prudential suggestions, as previously described. Postoperative Measures: The patient, the day after surgery, must be visited to check the general and local state of health (extension of edema, hematomas, dehiscence, any paresthesia and / or residual anesthesia). The sutures are removed six / ten days after surgery, after controlling the residual edema and gently disinfecting the wound. All the suture material must be carefully removed because any residue may have deleterious effect on the healing prognosis. Post-Rehabilitation Treatments: It is

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Istruzioni d'Uso
 Instruction of use
 Rev.2 del 09.11.2019

essential to establish regular and periodic monitoring of peri-implant tissues. This also includes the careful peri-implant bone assessment through clinical and radiological exams. Oral hygiene practiced by the patients must be monitored and encouraged. The intraosseous, transmucosal and prosthetic implant components must be of high quality. Over time, clinicians should monitor their effectiveness and efficiency. It is important to maintain the biological and bacterial balance in the various and different types of implant-prosthetic rehabilitations through an accurate patient home care. Mechanical tools, used to scaling and root planing, should be used with enough skill to avoid damaging at the transmucosal components. All informations about follow-up and control Rx, able to reveal the presence of disorders at the implant-prosthesis complex, should be recorded on the medical report. Packagin: the packaging, if kept intact, isolates the inner material from pathogens and preserve the sterile condition until the expiry date reported on the package. We accept no liability for the use of re-sterilizes implant. Caution: The expiry date reported on external and internal labels, refers to the sealed container properly stored. Do not use the product after the expiry date reported on the package. Disposal of the product must take place in accordance with the disposal protocols required by law as a special waste. Risk for Reuse of an Implanted Medical Device: Surface damage - Trigger cracks - Structural Collapse - Chronic inflammation - Surfaces wear - Localized corrosion of the implant – Toxicity. Medical Device to be used on adult patients only.

Manufactured by: AXELMED S.R.L. / Via della Liberazione, 58 / 20098 S. Giuliano Mil.se (MI) / Italy.

Diametro Impianto Implant Diameter	Connessione protesica Prosthetic Connection	Lunghezza (mm) Length (mm)
Ø 3,4 mm	Esagono interno Internal hex Ø 2,4 mm	9 - 11 - 13 - 15
Ø 3,8 mm		9 - 11 - 13 - 15
Ø 4,3 mm		7 - 9 - 11 - 13 - 15
Ø 5,0 mm		7 - 9 - 11 - 13 - 15
Ø 5,6 mm		7 - 9 - 11

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